

Extend™ system for Leksell Gamma Knife® Perfexion™

Article No : 1009633



Extend™ system enables fractionated treatments with Leksell Gamma Knife® Perfexion™. The non-invasive repeat fixation device is based on a vacuum assisted mouthpiece. Extend™ system is a fully integrated solution for Leksell Gamma Knife® Perfexion™.

Consisting of

- Extend™ Frame
- Extend™ Front Piece
- Patient Control Unit
- Extend™ License for Leksell Gamma Plan
- Extend™ PPS Hardware kit
- Reposition Check Tool
- RCT QA Tool
- Extend™ Vacuum Cushion 1L (5 pcs)
- Extend™ Vacuum Cushion 1.5L (3 pcs)
- Extend™ Mouth Piece kit
- Extend™ Patient Preparation kit
- Extend™ Dental kit
- Extend™ CT Indicator
- Extend™ CT Adapter Interface
- Extend™ Instruction For Use
- Extend™ Emergency Procedures
- Installation
- Training