


## Ear Plug Holders - One Pair

Article No : 1004312

A vertical strip on the left side of the page showing a close-up of the ear plug holders, which are small, clear, dome-shaped devices with three holes, attached to a metal frame.

Ear Plugs and Holders are optional devices that can be fitted to the side bars to facilitate alignment and temporarily support the frame on the patients head before it is secured in place. An ear plug fits into the holder through one of its three holes. **Note:** The Ear Plugs holder with article number 1004312 is compatible with the new Sterilization Tray for G frame with article number 1013578 (kit) and 1013823 (base tray), but is not compatible with the previous frame sterilization tray with article number 50151-03.