

Slotted Front Piece

Article No : 1002678



Slotted Front Piece can be used for 3-point fixation of Leksell® Coordinate Frame G enabling the surgeon to extend the access to the surgical area.

Included in the kit are 1 slotted front piece, 1 support, 1 locking piece for fixation post and 1 Screw 5x19 mm, Titanium.

By exchanging any of the standard front pieces for the Leksell Coordinate Frame G to the Slotted Front Piece, one of the two anterior Fixation Posts can be moved to any position closer to the mid-line allowing a total removal of the second anterior Fixation Post.