

Frame Fixation

Article No : 1000349



Frame Fixation offers a rigid and accurate fixation of Leksell® Coordinate Frame G to the operation table headrest. The three point fixation allows for high stability even if the patient is shaking.

Ensuring accuracy

The fixation design is the same as in our scanner adapters to ensure accuracy in treatment, i.e. the mechanical stress on the frame is the same during surgery as during imaging.

Frame Fixation is compatible with the Mayfield* or a similar headrest.

The interface between Leksell Coordinate Frame G and Frame Fixation is straightforward and easy to handle concerning attachment and fast release.

*Mayfield is a registered trademark of Ohio Medical Instrument Co, Inc