

Leksell Stereotactic System® Clamp

Article No : A1200-04



Leksell Stereotactic System® Clamp is used to attach the frame to the operation table headrest. The clamp provides an adjustable and rigid fixation to the operating table, facilitating the operative procedure and ensuring patient comfort. By using the standard screw, the Leksell Stereotactic System® Clamp is fixed to the Mayfield headrest's top joint.

After the patient is positioned on the table the clamp can hold any side of the frame and be adjusted to any desired position.

Mayfield is a registered trademark of Ohio Medical Instrument Co, Inc

Consisting of

- Frame Support, complete for LSS
- Knob, Leksell Stereotactic System® Clamp
- Allen Wrench 6 mm