

Reusable Fixation Screws, Al, Kit of 20 pairs

Article No : 1006581



Reusable Fixation Screw ensures quick and rigid skull fixation without compromising imaging quality. Reusable Fixation Screws are made of coated hard aluminum. The low density and non magnetic properties of the Reusable Fixation Screw makes it the optimum choice in scanning procedures such as MR and CT. Reusable Fixation Screws have female hexagonal sockets. This set consists of 20 pairs of reusable screws in various lengths. Expected life length is one year. The low density and non magnetic properties of the Reusable Fixation Screw mean that it is the optimum choice in scanning procedures such as MR and CT.

The screw has a female hexagon socket that allows for flush mounting with the fixation posts, allowing optimal treatment volume with Leksell Gamma Knife® and ease-of-use with tight MR head coils. Always use Insulated Fixation Posts during MRI scanning.

When mounting the frame with the reusable screws the instrument screwdriver is used (the instrument screwdriver is not included in the kit).

Reusable Fixation Screw can be sterilized using autoclaving, STERRAD and EtO methods.

Note: Every length can also be ordered separately in minimum order of screws.

Consisting of

- Reusable Fixation Screw, L= 20 mm
- Reusable Fixation Screw, L= 22,5 mm
- Reusable Fixation Screw, L= 25 mm
- Reusable Fixation Screw, L= 27,5 mm
- Reusable Fixation Screw, L= 30 mm
- Reusable Fixation Screw, L= 32,5 mm
- Reusable Fixation Screw, L= 35 mm
- Reusable Fixation Screw, L= 37,5 mm
- Reusable Fixation Screw, L= 40 mm
- Reusable Fixation Screw, L= 42,5 mm
- Reusable Fixation Screw, L= 45 mm
- Reusable Fixation Screw, L= 50 mm
- Reusable Fixation Screw, L= 55 mm
- Reusable Fixation Screw, L= 60 mm
- Reusable Fixation Screw, L= 65 mm
- Reusable Fixation Screw, L= 70 mm
- Reusable Fixation Screw, L= 75 mm
- Reusable Fixation Screw, L= 80 mm
- Reusable Fixation Screw, L= 90 mm
- Reusable Fixation Screw, L= 110 mm

- Reusable Fixation Screw, L= 110 mm
- Sterilization Case Fixation Screws

