

Fixation Screws, Titanium, Kit of 5 pairs

Article No : 811105



The screw has a hexagonal socket. This enables flush mount with the frame posts combined with maximum frame stability. The finest medical grade titanium enables extended lifetime.

Long lasting and durable fixation screws with a female hex design enables an expanded treatment volume in Leksell Gamma Knife® and increased ease-of-use in tight MR head coils. Always use Insulated Fixation Posts during MRI imaging.

The interaction between the screw and the post is improved thanks to the lubricating effect of the titanium.

When mounting the frame with the titanium screws the instrument screwdriver is used (the instrument screwdriver is not included in the kit).

Note: Every length can also be ordered separately in minimum order of two items.

Titanium Fixation Screw, Kit of 5 pairs

Consisting of

- Titanium Fixation Screw, 40 mm
- Titanium Fixation Screw, 45 mm
- Titanium Fixation Screw, 50 mm
- Titanium Fixation Screw, 55 mm
- Titanium Fixation Screw, 65 mm