

Insulated Fixation Posts Kit

Article No : 1006477



Insulated Fixation Posts effectively insulate the fixation screw from the post and provide a sterile, single-use interface between patient and frame. For MRI as well as CT scanning. Locking screws and locking pieces are included in the kit.

Insulation

Insulated Fixation Posts have sterile, single-use inserts for the fitting of the screws. The new design effectively insulates the fixation screw from the post, thereby minimizing the occurrence of radio frequency generated thermal effects when using MR scanners. Designed for 1.5 and 3.0 T MRI scanners.

Sterile Interface

These new posts are the result of the continuous research and development of the Leksell® Stereotactic Neurosurgery instrumentation in order to meet new needs at the forefront of stereotaxy. The risk of persistent crosscontamination is reduced when Insulated Fixation Posts are used together with Disposable Fixation Screws since the insert provides a single-use interface between the patient and the frame.

Increased Access

Insulated Fixation Posts have a rounded, slimmed top that both allows for more space in tight MRI head coils and increases the available treatable volume for Gamma Knife® surgery. A further improvement provided by the plastic inserts is a smoother interface between the metal post and the metal fixation screw. This increases the lifetime of the screws as well as ease-of-use when attaching the stereotactic frame.

Locking screws and locking pieces are included in the kit.

Consisting of

- Insulated Fixation Posts, Anterior - One pair
- Insulated Fixation Posts, Short Posterior - One pair
- Insulated Fixation Posts, Long Posterior - One pair
- Disposable Inserts