


## MR-Adapter for Siemens CP Head Array

Article No : 307031



Designed for Siemens CP Head Array coils. Inner Diameter 268 mm. Each adapter comprises a base which attaches to the MR scanner table and provides a rigid platform for the frame support. The MR adapter is constructed from a nonmagnetic material. The exact appearance of the adapter and the means by which it is attached to the scanner table vary according to the type of scanner. A range of MR adapters to suit different designs of MR scanners is available from Elekta.

For further details on which adapter you require for your scanner, please contact Elekta.