

MR-Adapter for Siemens Conical Head-Coil

Article No : A0420-24



This MR adapter is designed for Siemens Magnetom with circular polarized conical head-coil, shell type opening. Inner diameter of 268 mm. Each adapter comprises a base which attaches to the MR scanner table and provides a rigid platform for the frame support. The MR adapter is constructed from a nonmagnetic material. The exact appearance of the adapter and the means by which it is attached to the scanner table vary according to the type of scanner. A range of MR adapters to suit different designs of MR scanners is available from Elekta.

For further details on which adapter you require for your scanner, please contact Elekta.