

MR-Adapter Universal (Short)

Article No : A0420-29



This MR adapter can be an option for certain MR scanners when no other dedicated MR adapter is suitable. This is a universal MR adapter of a short design.

It can also be used as an option on certain X-Ray tables when no other suitable table fixation is available. Each adapter comprises a base which attaches to the MR scanner table and provides a rigid platform for the frame support. The MR adapter is constructed from a nonmagnetic material. The exact appearance of the adapter and the means by which it is attached to the scanner table vary according to the type of scanner. A range of MR adapters to suit different designs of MR scanners is available from Elekta.

For further details on which adapter you require for your scanner, please contact Elekta.