

CT-Table fixation for GE

Article No : 912965



This CT table fixation is designed for GE scanners with metal-free cradle and plastic latch for headrest attachment, eg. Pro-Speed, SX, HiSpeed, LightSpeed, Lemage SX and GE Revolution. The CT table fixation is fitted to the scanner table, usually in place of the headrest. It is designed to hold the CT adapter, coordinate frame, CT indicator and the skull of the patient firmly in the correct position and orientation within the scanner.

A range of CT table fixations to suit different designs of CT scanners is available from Elekta. This table fixation is designed for GE scanners with metal-free tables.

For further details on which table fixation you require for your scanner, please contact Elekta. This product is the same as product 1040043 registered under the Medical Device Regulation.