


CT-Table fixation for Philips

Article No : 1002031



This table fixation is designed for Philips Brilliance CT and IQon Spectral CT using a smooth top slotted for headrest fixation. The CT table fixation is fitted to the scanner table, usually in place of the headrest. It is designed to hold the CT adapter, coordinate frame, CT indicator and the skull of the patient firmly in the correct position and orientation within the scanner.

A range of CT table fixations to suit different designs of CT scanners is available from Elekta.

For further details on which table fixation you require for your scanner, please contact Elekta. This product is the same as product 1040040 registered under the Medical Device Regulation.