

Table Fixation for X-Ray, Flat 210 mm

Article No : A0440-02



Designed for flat table head end (GE DF 5000, Hitachi DH 158HM), width 210 mm. Table Fixation for X-Ray, Flat 210 mm The X-Ray Table Fixation is fitted to the table of the X-Ray machine. The table fixation is designed to hold the coordinate frame, indicator and the skull of the patient firmly and to ensure that position and orientation can be reproduced in subsequent scans.

A range of X-Ray table fixations to suit different designs of X-Ray tables are available from Elekta. This X-Ray Table Fixation is designed for flat table head end (GE DF 5000, Hitachi DH 158HM), width 210 mm.

For further details on which table fixation you require for your table, please contact Elekta.

X-Ray machines manufactured and/or supplied by GE Medical Systems or Hitachi may be trademarked and otherwise protected by GE Medical Systems or Hitachi. Elekta does not manufacture and distribute and can therefore assumes no liability whatsoever for the performance of this product.