

Table Fixation for X-Ray, Siemens Curved

Article No : A0440-01



Designed for curved table head end (Siemens), width 226 mm. Table Fixation for X-Ray, Siemens Curved The X-Ray Table Fixation is fitted to the table of the X-Ray machine. The table fixation is designed to hold the coordinate frame, indicator and the skull of the patient firmly and to ensure that position and orientation can be reproduced in subsequent scans.

A range of X-Ray table fixations to suit different designs of X-Ray tables are available from Elekta. This X-Ray Table Fixation is designed for curved table head end (Siemens), width 226 mm.

For further details on which table fixation you require for your table, please contact Elekta.

X-Ray machines manufactured and/or supplied by Siemens Medical Systems may be trademarked and otherwise protected by Siemens Medical Systems. Elekta does not manufacture and distribute and can therefore assume no liability whatsoever for the performance of this product.