

Disposable Biopsy Needle kits, set of 6

Article No : 911933



These needles are an alternative to reusable biopsy needles addressing concerns regarding cross-contamination. The kit includes a disposable side-cutting biopsy needle with an 8.0 mm needle window and an aspiration tube. The kit also includes a disposable Stop and Guide Insert to further reduce the risk of cross-contamination.

Product features

- Sterile
- Disposable
- Working length of 190 mm
- 8 mm side cutting window
- Inner and outer cannula with polished, rounded closed ends
- 14 gauge nominal diameter (2.11 mm)
- Including aspiration tube, disposable stop and guide
- Luer lock hubs provide a highly efficient seal during aspiration

The kit includes:

6 Biopsy needles with inner and outer cannulas in stainless steel
6 Biopsy needle guides
6 Biopsy needle stops
6 Aspiration tubes