

Universal Guide tube 10/13 mm before target point

Article No : MER230880



Universal Guide tube for MicroMacro or Micro Electrode and Medtronic DBS lead. The guide tube is fixated in the MicroGun and is not moving with the electrode when advanced towards the target. The inner guide tube for the MicroMacro or Micro Electrodes ends 10 mm before target and the outer guide tube for the Medtronic DBS lead ends 13 mm before target. Inner guide tube Outer Diameter = 1.40 mm, Inner Diameter = 1.0 mm, Length = 206 mm. DBS guide tube Outer Diameter = 1.65 mm, Inner Diameter = 1.45 mm, Length = 203 mm. Packed sterile, single use. Set of 5, individually packed.

Only for DBS electrodes longer than 300 mm.

Manufactured by inomed (230880).