


## Leksell Stereotactic System® Planned Maintenance on demand (Sweden)

Article No : 914385



Maintain the functionality and ensure the accuracy and overhauling of Leksell Stereotactic System to a good working condition. Functionality and performance are verified and certified. The parts are sent to Elekta in Sweden for service. Leksell Stereotactic System® Planned Maintenance on Demand

Provides the customer with a service of Leksell Stereotactic System® (LSS) where functionality and performance are verified and certified. Planned Maintenance can be executed once every year. Planned Maintenance is performed on contracted Leksell Stereotactic System® to the applicable extent. The service will be documented with a certificate. If necessary a service report will be provided describing outstanding service needs and or recommended parts replacement to sustain equipment performance.

Planned Maintenance Sweden includes: Cleaning, inspection of system, lubrication, refilling of engravings if needed, accuracy check and replacement of front piece screw, lock screw on clamp, nylon stop screw, knob and lock screw for X-axis.

Leksell Stereotactic System® parts covered by this service are:

- Leksell® Coordinate Frame G
- Leksell® Multi Purpose Stereotactic Arc including Arc supports with slides

For further information regarding this service please contact your local Elekta representative