

Radiation Phantom with Cassettes

Article No : A0202-02



The Radiation Phantom with Cassettes is used for calibrating the absorbed dose rate of Leksell Gamma Knife®. Radiation Phantom with Cassettes The phantom is divided into two halves between which a cassette can be inserted. The cassettes allow film dosimetry, TLD and dose measurement by means of a small ionization chamber. The cassettes are made of polystyrene and are 10 mm thick.

The material is selected to respond within 1% to ^{60}Co radiation as water would regarding absorbed dose.