


Leksell GammaPlan® DICOM-RT

Article No : 1016994



DICOM-RT provides an exclusive possibility to share treatment information with any software supporting DICOM-RT (Treatment Planning Systems, PACS).

It includes:

IMPORT of DICOM-RT STRUCT allows displaying any user-defined region/volume created on a DICOM-RT system. Volumes can be target volumes, organ at risks, isodoses or any other region of interest defined by the user.

EXPORT of treatment data via DICOM-RT STRUCT and DICOM-RT DOSE allows sharing information with any compatible systems to perform dose comparison or dose addition.