

## Leksell Stereotactic System® Planned Maintenance on demand (local)

Article No : 1000609



Maintain the functionality and ensure the accuracy and overhauling of Leksell Stereotactic System® to a good working condition. Functionality and performance are verified and certified.

Provides the customer with a service of Leksell Stereotactic System® (LSS) where functionality and performance are verified and certified. The service will be documented with a certificate. If necessary a service report will be provided describing outstanding service needs and or recommended parts replacement to sustain equipment performance.

Planned Maintenance includes: Cleaning, inspection of system, lubrication, accuracy check and replacement of front piece screw, lock screw on collar, nylon stop screw, knob and lock screw for X-axis.

Leksell Stereotactic System® parts covered by this service are:

- Leksell® Coordinate Frame G
- Leksell® Multi Purpose Stereotactic Arc including Arc supports with slides

For further information regarding this service please contact your local Elekta representative