

## Luer Connection Needle

Article No : 1020177



The Luer Connection Needle is part of the Injection/Aspiration Needle Kit. The needle is designed for use with Leksell Stereotactic System® for

diagnostic and therapeutic punctures, aspirations and evacuations.

The outer diameter of the needle is 1.5 mm.

The Luer Connection Needle contains a cannula with Luer lock and a stylet.

The working length of the needle is 190 mm.

This needle is replacing the previous product called Cannula with 3-Way stopcock and stylet 1.5 mm with part no 50067-01.

Manufacturer: Elekta Instrument AB