

## Elekta Stereotactic Drape

Article No : 1514116



Clear see-through and thin drape for fast, safe and sterile draping of Leksell Stereotactic System® and Leksell® Vantage™ Stereotactic System. The tailor-made design ensures a large and accessible sterile field. Equipped with an wide adhesive for fixation to the top of the patient's skull, which makes it optimal for a bilateral procedure. Designed with a drain pouch and instrument pouches as well as an adhesive on the anterior side of the drape to clear it above the patient face. There are two side slits for the arc supports for Leksell Stereotactic System. Set of 10 individually sterile drapes.

Legal manufacturer: Microtek medical, Inc. Distributed by: Elekta Instrument AB