

FirmFix™ 2

Article No : 1519663



The single use FirmFix™ pin attaches the Leksell Vantage™ Head Frame to the patient's head. The FirmFix™ are delivered in a box containing 20 sterile blisters with 2 FirmFix™ in each pouch. There are five sizes available, where FirmFix™ 1 has the shortest length and FirmFix™ 5 the longest length. The fixation pin is designed using a PEEK polymer and a Aluminium tip.

This product is the same as product 1040032 registered under the Medical Device Regulation.

This product is the same as product 1040032 registered under the Medical Device Regulation.