

Leksell Vantage Arc System

Article No : 1053958



Leksell® Vantage™ Arc System is designed to allow the neurosurgeon to select any desired trajectory to the predetermined brain target. Utilizing the proven Leksell Coordinate System and Center-of-Arc principle in an innovative click-on design and easy sterile draping together with Leksell® Vantage™ Head Frame. The arc assembly can be attached in a lateral-right or lateral-left position.

The arc system's purpose is to set planned target coordinates, trajectory and guide the instruments used during neurosurgery. It is used for diagnostic and therapeutic procedures, such as functional interventions and minimal invasive neurosurgical procedures within its intended use.