

MRI Head Coil Adapter, lower stop

Article No : 1054773



There are four different coil adapters which supports the patient's head within the frame inside the MR Head coils and prevent movement. The coil adapter is attached to the back of the frame.

The coil adapter's lower stop prevents the coil adapter from sliding too far into the head coil. and fits well inside a Siemens TIM Matrix Head coil

All the coil adapters have four rubber feet in the corners. This is to minimize the risk for movements when the patient has the frame, coil adapter and the fiducial box attached on the head during the MR imaging processing in the head coil.