

MRI Head Coil Adapter, upper stop

Article No : 1098534



MRI Head Coil Adapter, Upper Stop is used together with the Leksell Vantage Stereotactic System to support the patient's head inside a MR Head coil. The coil adapter attaches to the posterior part of the frame. This MRI head coil adapter is compatible with the following MR head coils:

Philips dStream T/R Head coil
Philips Transmit Receive Head Coil

This product replaces 1054774.

This MRI head coil adapter is tailored for the Philips 3T/1.5T coils. It is similar but slightly narrower, and thicker than the 'Universal' adapter. In addition, it has two 'towers' extending in the Z-direction that act as mechanical stops to prevent the MRI fiducial box from getting in contact with the inner wall of the coil.

This product is the same as product 1040075 registered under the Medical Device Regulation.