

MRI Head Coil Adapter, large

Article No : 1098530



MRI Head Coil Adapter, Large is used together with the Leksell Vantage Stereotactic System to support the patient's head inside a MR Head coil. The coil adapter attaches to the posterior part of the frame. This MRI head coil adapter is compatible with the following MR head coil(s):

GE Quadrature Head coil
GE 450/750 and HD split head coil
Siemens TxRx CP Head coil
Siemens CP Head Array Coil

This coil is designed to be used in any coil with a diameter larger than app 250 mm. It is similar to MRI Head Coil Adapter, Universal but with a larger diameter. Its rubber "feet" are also more stretched apart to optimize the stability of the patients head.

This product replaces 1054775.

This product is the same as product 1040080 registered under the Medical Device Regulation.