

Article No : TRT 0361

## ABC Full System



The Active Breathing Coordinator provides all features required to implement a breathholding technique within the clinical environment. This system is stand alone and can be easily utilized with a range of linear accelerators and CT scanners.

### Consisting of

- Control Computer with Control Software installed
- Flat Panel Patient Monitor
- CYBEX Transmitter/Receiver
- Transducer (pick-up assembly)
- Transducer (turbine)
- Balloon Valve
- Patient Control Switch
- Mirror & Mouthpiece Support System
- Mouthpiece & Filter one-time use kit (set of 40)
- ABC Cart System
- CAT 5 LAN Cable (please state required length in your order)