

Article No : TRT 3701

ABC R2.0 Full Upgrade Kit



This upgrade allows customers installed with Active Breathing Coordinator R1.02 to obtain the features included in R2.0.

Hardware

- Modified connection port on the mirror support system supports the Transducer turbine and pickup assembly and minimizes strain on the transducer cable.
- Color coding kit allows the user to identify the component parts of the respiratory system, to aid easy identification during assembly/ disassembly.
- Higher specification laptop with Windows XP, offering a more stable platform and ease of data.

Software

- Increased response time of the balloon valve, thereby improving the threshold accuracy capture.
- New user front end to allow for easier navigation and selection of control parameters.
- New record and pause functionality allowing specific selection of data captured in each session.
- New tools to allow retrospective analysis including:
 - View current session traces
 - Import up to 5 previous session traces
 - Align respiratory traces to specific breath-hold's
 - Manipulation of view, including zoom, pan and reset tools.
 - Volume and time measurement tools
- New format for Session History and Output File viewing.

Consisting of

- Mirror Support Connector with Strain Relief
- Laptop Computer with CD R/W Drive
- Coupler with Color Coding
- Release 2.0 Firmware
- Identification Color Coding Kit
- ABC Release 2.0 Software