

Article No : P10106-301

## Fraxion Thermoplastic Mask Nose (10 pcs)



The Fraxion Thermoplastic Masks provide a large number of immobilization points over the patient's head. The thermoplastic masks are available in nose (combination with mouthpiece) and head design. Both options are double-layered and in micro perforation format. It will be stored together with all patient-specific parts in the storage box for the whole treatment time.

The Fraxion Thermoplastic Mask Head and Nose is only compatible to the Fraxion Head Support Module ST (P10106-111) and the previous Fraxion Frame (P10106-100).