

Article No : P10102-753

BodyFIX 14 Stereotactic Localization Set

The BodyFIX 14 Stereotactic Localization Set includes the localizer to localize the isocenter of the tumor and the target positioner to position the patient correctly for the treatment. This set also includes a marking set and the localizer syringe.

Consisting of

- BodyFIX Localizer Syringe Set
- BodyFIX 14 Localizer STANDARD
- BodyFIX 14 Target Positioner STANDARD
- BodyFIX 14 Marking Set STANDARD

¹ All "BlueBAG BodyFIX 14" are indexed with Indexing lath (Standard).