

Article No : P10102-754

BodyFIX Basic Accessories Set

The BodyFIX Basic Accessories Set includes a stabilizing cushion 800x600 mm (31x24 Inch), a cushion OP portal 800x600 mm (31x24 Inch), a comfort cushion 500x200 mm (20x8 Inch) for the head and a comfort cushion 500x400 mm (20x16 Inch) for the upper part of the body.

Consisting of

- BodyFix Comfort Module 500x200 mm
- BodyFix Comfort Module 500x400 mm
- BodyFix Stabilizing Cushion 800x600 mm
- BodyFix Stabilizing Cushion OP-PORTAL

¹ All "BlueBAG BodyFIX 14" are indexed with Indexing lath (Standard).