

Article No : P10102-774

BodyFIX 14 TOTAL BODY Fixation Set WIDE

This set includes all necessary components to enable patient positioning especially for total body treatment.

Consisting of

- BodyFix Manifold Cushion 600 mm
- BodyFix Manifold Cushion 800 mm
- BodyFix Manifold Tube 2000 mm
- BodyFIX Manifold Connector STANDARD
- BlueBAG BodyFIX 14 Rectangular 900x2325 mm/117L

¹ All "BlueBAG BodyFIX 14" are indexed with Indexing lath (Standard).