

Article No : P10102-780

BodyFIX® 14 System Stereotactic Frame-based for Non-indexed Tabletops

This is a system for non-indexed tabletops. It's ideal for immobilization and accurate repositioning of the patient from planning to treatment delivery, using a stereotactic target positioner and localizer.

Consisting of

- BodyFIX 14 Base Plate Stereotactic STANDARD (2000 x 530 x 22 mm)
- BodyFIX 14 Marking Set STANDARD
- BodyFIX Cover Sheet TOTAL BODY 2700x1400 mm (set of 5)
- BodyFIX Hygienic Drape Manifold Tube (2x15000 mm)
- BodyFIX Adhesive Tape Set of 3
- BodyFIX Independent Indexing LINAC Couch Adapter (Set of 2)
- BodyFIX Independent Indexing CT Couch Adapter (Set of 2)
- BodyFIX Indexing Adapter P2 (set of 3)
- BodyFIX 14 Stereotactic Localization Set
- BodyFIX 14 TOTAL BODY Fixation Set
- BodyFIX 14 TOTAL BODY Fixation Set WIDE
- iBEAM Indexing Bars (Set of 3)
- BodyFIX Vacuum Pump P3

¹ All "BlueBAG BodyFIX 14" are indexed with Indexing lath (Standard).